

# SR-AM

## Head Holding Device (for Mice)

### Converting your stereotaxic instrument to a head holding device for mice.

The SR-AM head holding device can be attached to the base plate of SR series. Your stereotaxic instrument can be converted to an instrument for mice.



## Specification

<b>Accessories included</b>	<u>EB-3B</u> Ear Bars for Mice (in a pair) <u>EB-5N</u> Auxiliary Ear Bar for Mice Hex Wrench Mounting Knob
<b>Dimensions/Weight</b>	W300 × D190 × H88mm, 1.2kg

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2019