


RA-6N

Head Holding Adaptor for Rats

This adaptor is used to convert a mouse fixing clamp to a rat fixing clamp. The dovetail structure permits simple adjustment of the mouth and nose clamp position; the procedure can be smoothly and easily performed with one hand. This is also true of the clamping action using an auxiliary ear bar, which is easier than with conventional models.



- * An ear bar for mice ([EB-3B](#)) can be attached to an auxiliary ear bar for rats ([EB-4N](#)).
- * Without the auxiliary ear bar, an ear bar for mice ([EB-3A](#)) is required.

Specification

Accessory included	EB-4N Auxiliary Ear Bar for Rats
Dimensions/Weight	W20 × D85 × H55 ~ 75mm, 142g

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2019