

made to measure

TMR-02D

Timer



Features:

- ⇒ Generates square pulses of fixed amplitude (TTL, +5 V) and variable length
- ⇒ Adjustable PULSE length and DELAY
- ⇒ Times set by 5-digit thumb wheel switch
- ⇒ PULSE duration: 0.01 ms...999.99 ms
DELAY time: 0.01 ms...999.99 ms
RESOLUTION: 0.01 ms
- ⇒ Continuous cycling or preset number of cycles
- ⇒ Cycle counter with digital display
- ⇒ STOP and RESET function
- ⇒ START manually with toggle switch or remotely with TRIGGER signal (TTL, +5 V)
- ⇒ TRIGGER OUT signal (TTL, +5 V)
- ⇒ External power supply

For more information contact:

General:

npi electronic GmbH
Phone: +49-7141-9730230
Fax: +49-7141-9730240
sales@npielectronic.com
www.npielectronic.com

North America:

ALA Scientific Instruments
Phone: +1-631-393-6401
Fax: +1-631-393-6407
sales@alascience.com
www.alascience.com

Great Britain:

Scientifica Limited
Phone: +44-1825-749933
Fax: +44-1825-749934
info@scientifica.uk.com
www.scientifica.uk.com

Switzerland:

Science Products Trading AG
Phone: +41-43-4880561
Fax: +41-43-4880562
info@science-products.com
www.science-products.ch