

## GC-1/1.2/1.5 Glass Capillary

## Box of 200 (5 each to a vacuum packed bag).

These glass capillaries are pre-washed in an ultrasonic washing machine, and sterilized with ultraviolet ray. They are vacuum packed to protect against dust and dirt.





## GC-1/1.2/1.5 Glass Capillary

## **Specification**

Туре	<u>G-1</u>	<u>G-1.2</u>	<u>G-1.5</u>	<u>G-2</u>	<u>G-3</u>	<u>GD-1</u>	GD-1.2	<u>GD-1.5</u>	<u>GD-2</u>
Outer diameter (mm)	1	1.2	1.5	2	3	1	1.2	1.5	2
Inner diameter (approx. mm)	0.6	0.75	0.9	0.8	1.8	0.6	0.75	0.9	1.2
Length (mm)	90*	90*	90*	90*	90	90*	90*	90*	90*
Contents (pcs)	500	500	500	200	200	500	500	500	200
Material	Borosilicate								

Туре	GC-1	GC-1.2	GC-1.5	GDC-1	GDC-1.2	GDC-1.5	<u>G-100</u>	<u>G-100L</u>	<u>G-1000</u>	
Outer diameter (mm)	1	1.2	1.5	1	1.2	1.5	1	1	1	
Inner diameter (approx. mm)	0.6	0.75	0.9	0.6	0.75	0.9	0.75	0.75	-	
Length (mm)	90	90	90	90	90	90	90*	100*	90*	
Contents (pcs)	200	200	200	200	200	200	500	500	200	
Material	Borosilicate									

 $<sup>^{</sup>st}$  A custom length between 70mm and 180mm is also available by special order.

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2019